CA[®] POWER GRIP Set

**Bedienungsanleitung
Operation Manual
Mode d'emploi
Manual de Instrucciones**

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Pablo Echarri, D.D.S.

The CA® POWER GRIPS were designed by Dr. Pablo Echarri, with the help of Julio Blasco drawing the first drafts and José Fernández making the first real prototypes, both of them working in Ladent, Barcelona. The next step in developing the products were made in SCHEU-DENTAL, specially by Christian Scheu and Peter Stückrad, that have worked intensively to make it happen.

The CA® POWER GRIPS were designed to improve movements with the CA® CLEAR ALIGNER technique in a more effective and rapid way, reducing the steps of treatment.

The main concept in the MFM is to apply traction forces (elastics) in extrusion and rotation treatments at home and during the nights, and combining these forces with CA® CLEAR ALIGNER & attachments treatment during the day.

The helix shape design of the MDI and DMI allow the CA® CLEAR ALIGNER to control the mesio-distal inclination of the teeth offering more retention of the appliance over the teeth achieving the movements more easily.

The FALA and FALI increase the pressure of the CA® CLEAR ALIGNER in the gingival area of the teeth making possible the torque correction.

The BUB and BUM increase the vertical dimension minimizing the occlusal interferences and improving the movements of the teeth in cross-bite position.

The forms BUAL and BUAS are very useful to increase the torque during the incisors intrusion with more effective forces.

They can also maintain the intrusion forces while the patient is not using the CA® CLEAR ALIGNER, permitting molar extrusion.

All the CA® POWER GRIPS requires a special protocol of lab set-up.

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CA® POWER GRIP Set

REF 5338.1

Content: 18 silicone forms (9 versions, 2 pcs. each); including support

CA® POWER GRIP forms enable fabrication of individual attachments for use in the CA® CLEAR ALIGNER technique.

Set contains nine different forms which can be filled with light curing composite to be directly applied on the prepared tooth.

Curing with LED lamp.

(We recommend using Light Bond - Sealant Resin [02-0033] and Light Bond Resin - Paste [02-0031] or Tetric EvoFlow®) www.smile-dental.de

CA® POWER GRIP, 2 silicone forms each

1 MFM Multifunctional Forced Movement

For extrusion, mesial and distal rotation for use with CA® splints and elastics.

REF 5330.1

2 MDI Mesial / Distal Inclination

For setting upright frontal teeth, mesial /distal inclination.

REF 5331.1

3 DMI Distal / Mesial Inclination

For setting upright frontal teeth; distal / mesial inclination.

REF 5332.1

4 FALA Force Application LABIAL

For protrusion, retrusion and torque control up to 3°. For fabrication of activation surfaces on the teeth in the gingival area with vestibular positioning - for teeth 11, 13, 21, 23, 33, 43.

REF 5333.1

5 FALI Force Application LINGUAL

For protrusion, retrusion and torque control up to 3°. For fabrication of activation surfaces on the teeth in the gingival area with labial, lingual and palatal positioning, for teeth 12, 22, 31, 32, 41, 42.

REF 5334.1

6 BUB Build-Up Bicuspid

Round pad for premolars (3 x 3 mm) for occlusal bite elevation (approx. 2 mm).

REF 5335.1

7 BUM Build-Up Molar

Oval pad for molars (3 x 5 mm) for occlusal bite elevation (approx. 2 mm).

REF 5336.1

8 BUAS Build-Up Anterior Small

Occlusal bite plane for intrusion of upper frontal teeth (2.5 x 2 mm)

REF 5337.1

9 BUAL Build-Up Anterior Large

Occlusal bite plane for intrusion of upper frontal teeth (2.5 x 4 mm)

REF 5339.1

CA® POWER GRIP Support made of stainless steel, sterilisable

REF 5340.1

2 pcs.

Product Description:

The CA® POWER GRIP technique enables fabrication of individual attachments.

Nine different forms are available which can be filled with light curing composite and directly applied on the tooth which has been prepared with Sealant/Resin. Then proceed with light curing.

Application:

1. Preparation of teeth: Polish the tooth surface by means of a rotating brush or with polishing paste. Then proceed with rinsing and air-drying the teeth.
2. Etching: Apply etching gel or liquid on a mixed paper. Dry the teeth to be etched and apply the etching agent on the corresponding tooth surfaces, avoiding rubbing and instead dabbing it. This takes 30 seconds, in case of milk teeth 120 seconds. Avoid any contact of soft tissue and etching gel.
3. Rinsing and drying: Rinse teeth thoroughly with water spray, each tooth for 10 seconds. In case of using etching gel for 20 seconds in order to stop etching and to remove dissolved enamel particles.
4. CA® POWER GRIP Form: Place the required CA® POWER GRIP form on the supplied support.
5. Application of primer: Apply light-curing sealant/liquid on the etched enamel surface, e.g. Assure or Light Bond Sealant.
6. Application of paste: Fill the clean form completely with light-curing primer or e.g. Light Bond past medium or LCR paste. If necessary, use a small instrument in order to fill undercuts completely. Positioning the attachment: Position the form on the required place for the attachment by pressing it slightly.
7. Curing of paste: Curing is done within approx. 5 - 20 seconds, depending on lamp intensity and primer.
8. Removal of form: After material has completely cured, proceed with carefully removing the form moving it slightly laterally.

Cleansing / Disinfection:

- The used form is to be cleansed by customary mild disinfection solutions.
- Make sure the form is completely clean before disinfection. Do not use any disinfection solutions containing formaldehydes or high-concentrated bleaching agents, as these may lead to damage and discolorations. Please follow the instructions of the corresponding manufacturer.

To be used only by expert staff!

The above information is for the doctor's use only, not for the patient. It is the doctor who has to evaluate this information with regard to each special treatment case and to inform the patient regarding advantages, potential side effects and complications.

HINTS:

This information corresponds to our actual experience to the best of our knowledge. However, this information is of non-binding nature and it is subject to any modification within the context of technical progress or development. This information only describes the nature and function of our products and our services, providing no warranty claim. The buyer is held responsible for a careful examination of the products' functions and applications by qualified expert staff. The buyer has to follow existing laws and regulations on his own responsibility.

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SCHEU-DENTAL GmbH