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## Expansion set (4.0) for CA® photo support

The expansion set includes two fixing blocks for height adjustment of the camera for taking the following images:

1. Occlusal (upper or lower jaw)
2. Frontal (upper or lower jaw)
3. Lateral - left/right (upper or lower jaw)
4. Vestibular (upper and lower jaw)

**Note:** You require two fixing blocks for height adjustment for one jaw, i.e. two expansion sets will be required to take images of the upper and lower jaw!

The maximum height of stone-plaster models should not exceed 25 mm.

Labels (marked occlusal, frontal, left or right) should be attached to the fixing blocks for height adjustment - you'll find a model record sheet under the foam cover. You can print additional labels if required using the template at [www.scheu-dental.com/downloads/other](http://www.scheu-dental.com/downloads/other).

The expansion set also includes the CA® model holder with hexagon key. The dental stone models must be clamped in this model holder so that you can then take all the images mentioned above.

The accessories - CA® SMART Matrix, pins and sleeves - are required for transferring the CA® model holders into an articulator.

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## Instructions for using the expansion set (4.0) for CA<sup>®</sup> photo support

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### Adjusting the photo support:

Objective: Removal of the lower locking  
ring

1. Disconnect the USB device from the computer.
2. Initial situation: CA<sup>®</sup> photo support with locking ring.

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4. Carefully lift the complete camera holder vertically.

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5. Remove the locking ring (lower stop) as well.

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6. Note: Ball bearings are integrated in the camera holder on the guiding bar providing additional protection against twisting.

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7. Ensure that the ball bearings are correctly aligned with the groove when repositioning the camera holder. Tilting may dislodge the ball bearings.

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8. Attach the capture ring so that it is flush with the upper edge of the guide bar and fix it in place. This locked position later serves to fix the camera in place while the fixing block for height adjustment is being used.

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
9. Hook the elbow catch into the capture ring.

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
## Fitting the CA<sup>®</sup> model holder into an articulator

In order to prevent plaster penetrating the model holder during articulation, we recommend you fabricate a customised cover template made of DURAN<sup>®</sup> 0.75 mm (REF 3413). The CA<sup>®</sup> SMART Matrix (supplied with the expansion set) is reusable to facilitate production of the cover template.

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1. Place the CA<sup>®</sup> SMART Matrix as illustrated on the model platform - with the indentations facing upwards and not centrally!  
2. Pressurise DURAN<sup>®</sup> 0.75 mm - in case of corrections of the upper and lower jaw a cover plate must be fabricated for both jaws.

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
3. Mill out the three indentations with the CA<sup>®</sup> carbide bur.

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4. Place the side arms in an angled position up to the stop serving as basic positioning for articulation.

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5. Place the cover template on the model holder - the milled recesses must lie centrally over the blue sleeves.

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
6. Trim the base of the cover template with the SD foil scissors B to fit the respective articulator.

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7. Carefully insert the brass pins (REF 5313.850) into the blue sleeves (REF 5382.800) - ensure that the stop does not disappear in the sleeves!

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8. Note: If there is a reduction in the friction of the sleeves, the sleeves must be replaced. In doing so, ensure that the sleeves are inserted together with the pins into the hole. DO NOT glue in!

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9. Open the CA® model holder with the red hexagon key and clamp the models. Note: A longer screw (REF 5315.012) is available for bigger models.

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10. Check the space available between the upper cover template and articulator.

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11. Fit the upper model holder with dental stone in the articulator.

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12. After curing, turn the articulator (on its head) and position the clamped upper jaw model with check-bite.

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13. Fit the lower jaw model holder with applied cover template in the articulator.  
In illustration 13 + 14, the pattern mask has already been reduced according to the dimensions of the CA® SMART matrix.

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14. For all subsequent work steps with the CA® SMART 4.0, the models can be simply removed from the articulator and reset - the cover template must not be removed, otherwise it will not be possible to position the model correctly!

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
15. We recommend fixing the reduced pattern mask to the base plaster using CA® CYANO VENEER (REF 5302).

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Images in different views can now be taken as required using the CA® model holder in combination with the CA® SMART 4.0 Software.

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## Using the two fixing blocks for height adjustment in combination with the CA<sup>®</sup> model holder

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1. Place the fixing block for height adjustment on the rear edge of the base plate.

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2. To take occlusal images, push the fixing block for height adjustment to the required side.

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3. Lower the camera holder down to the stop plate.

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4. Place the CA<sup>®</sup> model holder with prepared dental stone model on the base plate.

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
5. Set the exact measurement level by turning the plate.

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6. Finally, secure the height position required by tightening the counter nuts and maintaining the position of the stop plate.

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
17. Set the exact height by turning the plate to the measurement level and then tightening the counter nuts.

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18. Proceed in the same way for the right side.

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
19. Note: If larger models are used and the height adjustment range is insufficient, plates and counter nuts can be screwed on upside down.

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## Using the CA<sup>®</sup> CHECKER in combination with the expansion set (4.0) for CA<sup>®</sup> photo support

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In order to precisely determine the measurement values (angulation and torque) with the CA<sup>®</sup> CHECKER before and after set-up, you must ensure that the CA<sup>®</sup> model holder of the expansion set lies flat on the model support. To do this, take the following steps:

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1. Cut off the heads of the blue sleeves flush with the metal base using a scalpel.

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2. Carry out a check on the model support.

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3. Take care to handle the CA<sup>®</sup> model holder on the model support correctly during measurement.

SCHEU-DENTAL GmbH  
[www.scheu-dental.com](http://www.scheu-dental.com)

phone +49 2374 9288-0  
fax +49 2374 9288-90

Am Burgberg 20  
58642 Iserlohn · Germany

DE/GB 1.000/05/16 G REF PM0141.01