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**How to fabricate a CA<sup>®</sup> CLEAR ALIGNER splint using  
VECTOR<sup>®</sup> 30 screws and the CA<sup>®</sup> CHROMIUM arch**

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# How to fabricate a CA<sup>®</sup> CLEAR ALIGNER splint using VECTOR<sup>®</sup> 30 screws and the CA<sup>®</sup> Arch CHROMIUM

Correction of malocclusion with a CA<sup>®</sup> CLEAR ALIGNER splint gets even more economic, effective and time-saving in comparison to the conventional CA<sup>®</sup> Set Up technique when using the specific VECTOR<sup>®</sup> 30 screw and the CA<sup>®</sup> Arch CHROMIUM.

Relevant literature:  
CA<sup>®</sup> Book Volume 2 Advanced  
(Chapter: VECTOR<sup>®</sup> 30 Screw with CA<sup>®</sup> Arch CHROMIUM  
[English: REF 5309 / German: REF 5307 / Spanish: REF 5306])

## Working Material:

REF 2484	VECTOR <sup>®</sup> 30
REF 5317	CA <sup>®</sup> CHROMIUM Arch
REF 3207	ISOFOLAN <sup>®</sup>
REF 3415	DURAN <sup>®</sup> 1.0 mm
REF 3259	DURASPLINT <sup>®</sup> monomer
REF 3258	DURASPLINT <sup>®</sup> polymer
REF 5302	CA <sup>®</sup> Cyano Veneer Fast
REF 3381	DIMO <sup>®</sup> PRO
REF 3226	BIOPERM-Trimmer

## Additionally recommended material:

REF 3442	SIL-KITT red
REF 8450	Ortho arch bender
REF 5300	CA <sup>®</sup> pipettes
REF 3220	Blocking out putty
REF 2041	CARAT <sup>®</sup> flat pliers
REF 5372	Set-Up manual saw
REF 5263	CA <sup>®</sup> grinding set
REF 3712	SD polishing brush
REF 5288	Diamond separating disc
REF 3714	SD high gloss polishing disc
REF 3710	SD polishing paste A
REF 3606	SD separating disc, serrated
REF 3460	SD foil scissors A
REF 3461	SD foil scissors B

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Cut off ISOFOLAN® at model rim and insert interdental cuts using a scalpel.

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Embed model into the pellets and thermoform DURAN® 1.0 mm. Then roughly cut it out using the SD foil scissors B.

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Push the two VECTOR® 30 screws onto the CA® arch CHROMIUM.

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For easier insertion the ends of the CA® arch CHROMIUM are conically tapered.

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
Reduce CA® arch CHROMIUM and bend the retention loops, making sure the retention ends can be used as selective supports on the splint in the molar region.

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If necessary, the lower part of the plastic space holder can be reduced, in order to make sure the screws are positioned in an right angle on the anterior teeth.

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Fix the CA® Arch CHROMIUM and the VECTOR® 30 screws with CA® Cyano Veneer Fast on the DURAN® material. Apply just one drop each to the opening of the bent retentions using the CA® pipette (see fig.)

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Cover selected areas on the splint using blocking-out putty for protection.

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Mix and apply DURASPLINT® according to the instructions and polymerize it in a pressure pot.

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Then roughly cut out the resin using the CA® carbide bur and mark the saw cuts with a pencil.

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Separate the material within the marked area using the set-up manual saw or separating disc. If necessary, leave the splint on the model during this process for stability reasons.

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Lift off the splint from the model and combine/extend the saw cuts between the screws with the HM tricutter fine (green).

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Shorten the splint using SD foil scissors A and finish it with BIOPERM trimmer.

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A dressed DIMO<sup>®</sup> PRO (see fig.) and SD high-gloss polishing mop in combination with SD polishing paste A are best suited for finishing the rims. Ensure to use them at a low speed (5.000 to 10.000 rpm).

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
Finally conventionally polish the areas where you applied the resin.

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Finished CA<sup>®</sup> CLEAR ALIGNER splint for the oblique segment expansion according to Dr. Pablo Echarri.

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**Important:**

Make sure not to exceed the maximal tooth movements in the CA® CLEAR ALIGNER technique acc. to Dr. P. Echarri (CA® Book volume 2, Advanced, pages 52/53).

**Hint:**

We recommend our course programm on this topic.

For further information and details on our courses please visit our website:

→ [www.scheu-dental.com](http://www.scheu-dental.com) → SCHEU-ACADEMY → Topics: CA® CLEAR ALIGNER

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REF 2484, REF 5317

**CA® CLEAR ALIGNER  
with VECTOR® 30 screws  
and CA® Arch CHROMIUM**

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REF 2487

**CA® CLEAR ALIGNER  
with VECTOR® 45 screw  
acc. to Prof. Hinz**

**TFA Schiene  
mit TFA Kit  
n. Dr. Echarri**

REF 5329

**TFA Splint  
with TFA Kit  
acc. to Dr. Echarri**

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**CA® CLEAR ALIGNER  
mit CA® Distalizer Kit  
n. Dr. Echarri**

REF 5319

**CA® CLEAR ALIGNER  
with CA® Distalizer Kit  
acc. to Dr. Echarri**