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How to fabricate a CA[®] CLEAR ALIGNER splint using the VECTOR[®] 45 / VECTOR[®] 45 M expansion screw acc. to Prof. Hinz

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How to fabricate a CA[®] CLEAR ALIGNER splint using a VECTOR[®] 45 / VECTOR[®] 45 M screw acc. to Prof. Hinz

By incorporating specific VECTOR[®] expansion screws in CA[®] CLEAR ALIGNER splints, treatment of malocclusions (e.g. crowdings or rotations of single teeth or groups of teeth) gets even quicker and more efficient and economical in comparison with the conventional CA[®] Set-up technique.

Thanks to double guiding bars, CA[®] CLEAR ALIGNER splints with an incorporated VECTOR[®] 45 / VECTOR[®] 45 M screw ensure secure and torsion-free fan-type expansion.

This manual describes **two different methods** of inserting the VECTOR[®] 45 / VECTOR[®] 45 M screw in a CA[®] CLEAR ALIGNER splint - the conventional process of polymerization in a **pressure pot** or alternatively **pressure moulding** which does not require using a pressure pot.

Note:

For the pressure moulding method it is mandatory to use a thermoforming machine with positive pressure in order to guarantee complete polymerization of the used resin (BIOCRYL[®] RESIN).

We recommend the following specialized literature for further information:

CA[®] Book Volume 2 Advanced

English: REF 5309 / German: REF 5307 / Spanish: REF 5318

Important:

Make sure not to exceed the maximal tooth movement in the CA[®] CLEAR ALIGNER technique acc. to Dr. Pablo Echarri (CA[®] Book volume 2, Advanced, pages 52/53).

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Working material for polymerization in a pressure pot:

REF 2487	VECTOR® 45
REF 2488	VECTOR® 45 M
REF 3207	ISOFOLAN®
REF 3407	CA® foil hard (0.75mm)
REF 8136	STEADY-RESIN monomer
REF 8134	STEADY-RESIN S polymer or
REF 8135	STEADY-RESIN M polymer or
REF 8140	STEADY-RESIN VARIO polymer
REF 3701	SD Cyano Veneer Fast
REF 3703	SD pipette

Additionally recommended material:

REF 3442	SIL-KITT red
REF 5288	diamond separating disc
REF 2041	CARAT® flat pliers
REF 5372	Set-up manual saw
REF 3606	SD separating disc, serrated
REF 3265	LC bur C
REF 3460	SD foil scissors A
REF 3461	SD foil scissors B
REF 3226	BIOPERM trimmer
REF 3382	DIMO® brown, coarse
REF 3381	DIMO® PRO
REF 3710	SD polishing paste A
REF 3711	SD polishing paste B
REF 3712	SD polishing brush
REF 3713	SD felt wheel
REF 3714	SD high gloss polishing disc
_____	Sticky wax (e.g. Supradent)

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Working material for pressure moulding:

REF 2487	VECTOR® 45
REF 2488	VECTOR® 45 M
REF 3207	ISOFOLAN®
REF 3407	CA® foil hard (0.75mm)
REF 3218	BIOCRYL®-RESIN monomer
REF 3221	BIOCRYL®-RESIN polymer

Additionally recommended material:

REF 3442	SIL-KITT red
REF 5288	diamond separating disc
REF 2041	CARAT® flat pliers
REF 5372	Set-up manual saw
REF 3606	SD separating disc, serrated
REF 3265	LC bur C
REF 3460	SD foil scissors A
REF 3461	SD foil scissors B
REF 3226	BIOPERM trimmer
REF 3382	DIMO® brown, coarse
REF 3381	DIMO® PRO
REF 3710	SD polishing paste A
REF 3711	SD polishing paste B
REF 3712	SD polishing brush
REF 3713	SD felt wheel
REF 3714	SD high gloss polishing disc
_____	Sticky wax (e.g. Supradent)

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
POLYMERIZATION IN A PRESSURE POT:

(Please refer to the pressure moulding method described on page 7)

Block out the model in the area of crowdings and undercuts with SIL-KITT and place it on the model platform of the pressure moulding device.

Thermoform ISOFOLAN® according to the manufacturer's instructions, reduce the foil at the model rim and make interdental incisions.

Embed the model into the pellets at about 4-6 mm below the cervical rim.

Pressurize a CA® foil hard, making sure the CA® logo is legible and placed on the last molars.

Now roughly cut out the mould while it is still on the model using the SD foil scissors B.

Reduce the lower part of the VECTOR® 45 /VECTOR® 45 M spacer with a diamond separating disc. **Note:** Make sure not to squeeze the spacer when cutting with scissors or reamer. If required, adjust the guiding pins with a reamer to the width of the jaw. **Note:** This will affect the expansion!

In order to fix the screw, apply a few drops of SD Cyano-Veneer Fast with the SD pipette to the desired position.

Note: The screw should be positioned behind the front beneath the teeth on the midline, making sure not to affect the counter bite!

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Mix and apply STEADY-RESIN polymer and monomer either in the „salt and pepper“ or „modelling“ technique to the area around the screw. Make sure the material transitions are smooth!
 Now polymerize the resin as usual in a pressure pot.

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Remove the plastic spacer with a pair of pliers.

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Carry out the separation cut to the splint while it is still on the model using the Set-up manual saw taking care not to damage the VECTOR® 45/VECTOR® 45 M screw.

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
Carefully lift off the splint from the model.

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Now roughly cut out the CA® CLEAR ALIGNER splint.

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
Then carry out a complete separating cut downwards and expose the adjustment unit of the screw with a scalpel beginning at the basis.

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
Smooth the resin areas using the DIMO® wheel (REF 3382) and a bur and clean the saw cuts.

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Finish the splint as usual. We recommend using the BIOPERM trimmer and the DIMO®PRO wheels.

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The materials listed on page 3 are ideally suited for high-gloss polishing.

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CA® CLEAR ALIGNER splint with visible activation arrow.

15 Fertige Oberkiefer CA® CLEAR ALIGNER Schiene.

Finished upper CA® CLEAR ALIGNER splint.

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
**POLYMERISATION BY PRESSURE
MOULDING:** Block out the model in the
area of crowdings and undercuts with
SIL-KITT. Mark the midline and drill a
hole for positioning the lower part of the
VECTOR® 45 /VECTOR® 45 M spacer,
taking care the screw is positioned as
close as possible to the model **and to
consider sufficient space for the lower
jaw front.**

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Place the model on the model
platform. Thermoform
ISOFOILAN® according to the
manufacturer's instructions,
reduce the foil at the model rim
and make interdental incisions.

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
Reduce and round off the upper part
of the VECTOR® 45/VECTOR® 45 M
spacer with a diamond separating
disc.

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Note: Make sure not to squeeze the
spacer when cutting with scissors or
cutter.

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
If required, adjust the guiding pins
with a cutter to the width of the
jaw.

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
Note: This will affect the expansion!

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Apply a few drops of sticky wax
to the drilling hole with the
ISOFOILAN® foil being on the
model, properly position the screw
within the hole and close to the
palate and fix it.

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The model should be positioned
with the front teeth facing the
pressure chamber. Then embed the
model into the pellets, making sure
the VECTOR® 45/VECTOR® 45 M
screw has as little undercuts as
possible. The model front should
be placed more deeply into the
model pot and the pellets should
be filled up only to the model rim.

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Important:

In order to avoid air inclusions, make sure to cover the palatal area with pellets as well.

The palatal area should be cut off up to approx. 1 cm before the positioned screw.

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Fix a CA® foil hard and mix BIOCRYL® RESIN until it has a creamy consistency. Apply a few drops of BIOCRYL® resin to the metal parts during the heating time of the CA® foil.

During the last 5 seconds of the heating time more resin should be added generously. Once the heating time has ended, close the pressure chamber and pressurize.

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The model should stay under pressure in the unit for at least 10 minutes. Ignore any acoustic signals during this period, making sure not to evacuate or to switch off the machine.

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Cut out material coarsely and separate it along the midline using a manual saw or a separating disc while the model is still on the model for stability reasons.

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
Expose the CA® foil in the area of the screw spacer with a bur.

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Lift off the splint from the model and remove the lower part of the spacer using a pair of pliers. Then reduce the rims to the required length and complete the separating cut.

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

The expansion screw can be pre-activated by a few turns so that the upper spacer can be easily removed.

Hint: Using the expansion key gets simpler if its end is bent a little.

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Finish and polish the splint as usual.

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We recommend using the material listed on page 3.

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

VECTOR® 45: ca. 3,0 mm

VECTOR® 45 M: ca. 6,0 mm

**Activation hint
for VECTOR® 45/VECTOR® 45 M
screw:**

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

→ www.scheu-dental.com → SCHEU-ACADEMY → Seminar → Themen: CA® CLEAR ALIGNER

Hint:

Please note this instruction does not replace attending our technical CA® courses.

We also recommend participating in the following additional course:

„Course No. 17 - module 5 - CA® CLEAR ALIGNER splints with specific VECTOR® expansion screws“.

For further information and details on the course please visit our website:

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More instructions on how to fabricate aligners and splints with different VECTOR® screws:

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REF 2484, REF 5317

**CA® CLEAR ALIGNER
with VECTOR® 30 screws
and CA® Arch CHROMIUM**

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REF 2487/2488

**CA® CLEAR ALIGNER
with VECTOR® 45/45 M screw
acc. to Prof. Hinz**

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REF 5329

**TFA Splint
with TFA Kit
acc. to Dr. Echarri**

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REF 5319

**CA® CLEAR ALIGNER
with CA® Distalizer Kit
acc. to Dr. Echarri**