IMPRIMO® LC Splint

Produktinformation

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IMPRIMO® LC Splint

Product information

Description:	IMPRIMO® LC Splint is a light-curing methacrylate-based resin that is specifically engineered for production of high-precision occlusal splints and for surgical drilling guides. Photo-polymerization is done by means of image projection systems (385 nm). IMPRIMO® LC Splint has been tested in terms of bio-compatibility and meets both the applicative and mechanical requirements. IMPRIMO® LC Splint should be used in building processes with a z-resolution of 25 - 100 microns. A minimum wall thickness of 1.5 mm is required for reasons of stability. Printer parameters available for Asiga MAX™ and Asiga PRO 4K.
Processing instructions:	Make sure IMPRIMO® LC Splint is tempered prior to use to 23 °C - 30 °C. Fill the material into the tray wearing your personal protective equipment (safety glasses, gloves and face mask). Small air bubbles can be removed carefully with a spatula. In case of contamination of the material such as dirt or fragments, the material can be filtrated and repeatedly mixed due to its low viscosity. Leave the material to stand for approx. one hour to avoid any air inclusions.
Contra-indications:	IMPRIMO® LC Splint should not be used for any purpose other than the generative CAD / CAM fabrication of high-precision occlusal splints and surgical drilling guides. IMPRIMO® LC Splint is not suitable for printing aligners. IMPRIMO® LC Splint is a medical device class IIa.
Post-processing:	Upon completion of the printing process, the objects shall be removed from the model platform. Uncured resin residues can be removed completely using Isopropanol and a soft brush or can be resolved in an appropriate cleaning device such as IMPRIMO® Clean or RS wash using IMPRIMO® Cleaning Liquid. Dry parts, until there are no residues of isopropyl alcohol (IPA) or equivalent cleaning agent.
Post curing:	Complete curing of the objects is done in the curing unit IMPRIMO® Cure or RS cure. The use of protective gas or vacuum prevent the creation of a dispersion layer is crucial. A wavelength of 385 nm is required. The polymerization process in the IMPRIMO® Cure takes 3 minutes in a nitrogen environment. In the RS cure, polymerisation takes 4 minutes in a vacuum. The corresponding parameters are already stored in the IMPRIMO® Cure and RS cure devices. Any deviation from the described post-curing process can lead to modifications of the mechanical properties or the material colour.
Shelf life / storage:	Store the material at room temperature (18 °C to 28 °C) and protect it from light. Any exposure to light might damage the material. Carefully seal the container immediately after use.
Hazard warnings:	May cause allergic skin reaction. May be harmful to aquatic organisms with long lasting effects. Attention: Polymerized resins are chemically resistant, avoid stains on clothes.
Disposal:	According to local authority prescriptions. Cured objects can be disposed of with other household waste. Un-polymerized material has to be disposed of at a waste collection point. All instructions on the processing of our materials - in oral, written or practical form - are given to the best of our knowledge and must be read as references. Material use and processing is beyond our control and lies within the responsibility of the user.
Technical data:	For professional use only. Density: approx. 1.1 g/L Viscosity: approx. 0.7 Pa s Flexural strength: 93 MPa Elongation at break: 10 % Solubility: 1.4 µg/mm ³ Water absorption: 24 µg/mm ³ Shore hardness: 80 D Comply: biocompatibility (irritation and delayed-type hypersensitivity; genotoxicity, carcinogenicity and reproductive toxicity; systemic toxicity; cytotoxicity) *These data result from measurements of a representative sample. Contains: Alkoxylated bisphenol-A dimethacrylate, initiators, stabilizers, dyes and pigments. Material Safety Data Sheet: www.scheu-dental.com/en/service

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SCHEU-DENTAL GmbH