IMPRIMO® LC Splint Produktinformation

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IMPRIMO® LC Splint Product Information

Description:

IMPRIMO® LC Splint is a light-curing resin at the base of methacrylate, specifically engineered for fabrication of high-precision dental occlusal splints and drilling templates. Use in printers with DLP technology using LED 385 nm. A minimum wall thickness of 1.5 mm must be kept. This means that the material is not suitable for printing aligners.

Processing instructions/ application:

Fill the material into the tray wearing your personal protective equipment (safety glasses, gloves and face mask). Small air bubbles can be removed carefully with a spatula. In case of contamination of the material such as dirt or fragments, the material can be filtrated and repeatedly mixed due to its low viscosity. Leave the material to stand for approx. one hour to avoid any air inclusions.

Post-processing:

Upon completion of the printing process, the objects have to be removed from the model platform. Un-cured resin residues can be completely removed using isopropyl alcohol and a soft brush.

Post exposure:

Final curing of the objects is done in an appropriate light oven such as for example the pro3dure medical CD-1. It is absolutely necessary to use protective gas in order to prevent the creation of a dispersion layer. A wavelength of 385 nm is required. The polymerization process in the light oven takes 4 minutes.

Any deviation from the described manufacturing process can lead to modifications of the the mechanical properties or the colour of the material. The biocompatibility of the material can only be guaranteed upon complete polymerisation. For oral use, the material must be completely cured.

Shelf life/storage:

Store the material at room temperature (18 °C to 28 °C) and protect it from light. Any exposure to light might damage the material. Carefully seal the container immediately after use.

Hazard warnings:

Prolonged skin contact with un-cured material and inhalation of monomer vapours can lead to irritations. In single cases, allergic reactions to certain ingredients of IMPRIMO® LC Splint might occur, for more details, please refer to our safety data sheets on IMPRIMO® LC Splint. In the event of eye contact, rinse thoroughly with water and seek medical advice. In the event of skin contact, rinse thoroughly with running water and soap.

Attention: Polymerized resins are chemically resistant, avoid stains on clothes.

Disposal:

According to local authority prescriptions. Cured objects can be disposed of with other household waste. Unpolymerized material has to be disposed of at a waste collection point.

All information on processing our material - verbal, written or practical - is given to the best of our knowledge and must be read as references without obligation. Material use and processing is beyond our control and lies within the responsibility of the user.

For professional use only.

Material Safety Data Sheet: www.scheu-dental.com/en/service

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