SCHEU GROUP

Tiefziehtechnik
Thermoforming Technique
Technique de Thermoformage
Tecnica de termomoldeo
Tecnica di termoformatura

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
BIOCRYL® C

Indication

Immediate dentures, orthodontic expansion plates, retainer plates, denture bases.

Tips and tricks:

Use of ISOFOLAN® required, wire elements and screws can be polymerised using BIOCRYL®-RESIN, take care while polishing (brittle)

Available in:

2.0 / 3.0 mm x 125 mm Ø
clear

2.0 mm x 125 mm Ø
pink, red, blue, yellow-
transparent

Technical data

Material:	PMMA
E-modulus:	3,300 MPa
Hardness:	Rockwell 195
Tensile strength:	70 MPa
Impact strength:	20 KJ/m ²

CE

BIOPLAST®

Indication

Duplicating moulds, sports mouth guard, bruxism splints, positioner, OSAMU-Retainer®.

Tips and tricks:

Can be melted together using IR heater or heat gun. Thermal polishing possible. Cleaning and perforation required. Use of BIOPLAST® insulating agent recommended.

Available in:

1.0 / 1.5 / 2.0 / 3.0 / 4.0 / 5.0 mm
x 125 mm Ø

Technical data

Material:	Ethylene-vinyl acetate
E-modulus:	15 MPa
Hardness:	approx. 78–82 Sh. A
Tensile strength:	18 MPa
Impact strength:	no fracture

CE 0044

BIOPLAST® bleach

Indication

Soft pliable rubber like material, best suited for fabrication of bleaching splints.

Tips and tricks:

Use of BIOPLAST® insulating agent recommended.

Available in:

1.0 mm x 125 mm Ø

Technical data

Material: Ethylene-vinyl acetate
E-modulus: 20 MPa
Hardness: approx. 92 Sh. A
Tensile strength: 20 MPa
Impact strength: no fracture

BIOPLAST® Color/Multicolor

Indication

Coloured mouth guard, positioner

Tips and tricks:

Can be melted together using IR heater or heat gun. Thermal polishing possible. Cleaning and perforation required. Use of BIOPLAST® insulating agent recommended.

Available in:

3.0 mm x 125 mm Ø
Various colour combinations

Technical data

Material: Ethylene-vinyl acetate
E-modulus: 20 MPa
Hardness: approx. 78-82 Sh. A
Tensile strength: 18 MPa
Impact strength: no fracture

BIOPLAST® XTREME

Indication

Individual sports mouth guard with integrated reinforcement in the anterior region.

Tips and tricks:

Additional reinforcement can be achieved by laminating DURAN® and BIOPLAST®. Cleaning and perforation required, use of BIOPLAST® Insulating agent recommended.

Available in:

4.0 / 5.0 mm x 125 mm Ø
red, blue, black and clear

Technical data

Material: Ethylene-vinyl acetate
E-modulus: 15/20 MPa
Hardness: 80/90 Sh. A
Tensile strength: 18/20 MPa
Impact strength: no fracture

COPYPLAST®

Indication

Duplication moulds, bleaching splints, bracket transfer masks, copings, moulds for temporary crowns and bridges, space maintainer

Tips and tricks:

Insulation against all kinds of resins, also suitable for the production of PMMA duplicates.

Available in:

0.5 / 0.75 / 1.0 / 1.5 / 2.0 mm
x 125 mm Ø

Technical data

Material: Polyethylene
E-modulus: 175 MPa
Hardness: 42 Sh. D
Tensile strength: 9 MPa
Impact strength: 101 KJ/m²

COPYPLAST® C

Indication

Long-term splints and retainers, orthodontic splints.

Tips and tricks:

Can be trimmed with scissors or scalpel, edges can be smoothed by means of DIMO® wheel or silicone polisher.

Available in:

1.0 mm x 125 mm Ø

Technical data

Material:	Polypropylen
E-modulus:	450 MPa
Hardness:	68 Sh. D
Tensile strength:	19 MPa
Impact strength:	650 KJ/m ²

C € 0044

DURAN®

Indication

Miniplast splints, occlusal splints, interim and fluorisation splints, surgical guides, anti-snoring appliances, reinforced anterior section of mouth guards, temporary crowns and bridges.

Tips and tricks:

Occlusal adjustments can be done using DURASPLINT® without bonding agent, use of ISOFOLAN® required.

Available in:

0.4 / 0.5 / 0.625 / 0.75 / 1.0 / 1.5 / 2.0 / 3.0 mm x 125 mm Ø

Technical data

Material:	PETG
E-modulus:	2,200 MPa
Hardness:	78 Sh. D
Tensile strength:	53 MPa
Impact strength:	10 KJ/m ²

C € 0044

DURAN®+

Indication

Miniplast splints, occlusal splints, interim and fluorisation splints, surgical guides, anti-snoring devices, temporary crowns and bridges.

Tips and tricks:

Occlusal adjustments can be done using DURASPLINT® without bonding agent. Comes with integrated insulation foil that is removed after thermofor-
ming, use of ISOFOLAN® not required.

Available in:

0.4 / 0.5 / 0.625 / 0.75 / 1.0 / 1.5 / 2.0 mm
x 125 mm Ø ; with integrated protection
foil

Technical data

Material:	PETG
E-modulus:	2,200 MPa
Hardness:	78 Sh. D
Tensile strength:	53 MPa
Impact strength:	10 KJ/m ²

CE 0044

DURAN®+ A2 pd

Indication

Tooth coloured temporaries, splints for planning and demonstration for aligner treatments, friction sleeves for the telescope technique.

Tips and tricks:

Occlusal adjustments can be done using DURASPLINT® without bonding agent. Comes with integrated insulation foil that is removed after thermofor-
ming, use of ISOFOLAN® not required.

Available in:

0,5 / 0,625 / 0,75 mm
x 125 mm Ø ; with integrated protection
foil

Technical data

Material:	PETG
E-modulus:	2,200 MPa
Hardness:	78 Sh. D
Tensile strength:	53 MPa
Impact strength:	10 KJ/m ²

CE 0044

DURAN®+ white pd

Indication

Cosmetic splints

Tips and tricks:

Occlusal adjustments can be done using DURASPLINT® without bonding agent, comes with integrated insulation foil that is removed after thermoforming, use of ISOFOLAN® not required.

Available in:

0,5 / 0,625 / 0,75 / 1,0 mm
x 125 mm Ø ; with integrated protection foil

Technical data

Material:	PETG
E-modulus:	2,200 MPa
Hardness:	78 Sh. D
Tensile strength:	53 MPa
Impact strength:	10 KJ/m ²

CE

DURASOFT® pd

Indication

Miniplast splints, bruxism splints, occlusal splints, intra-oral anti-snoring devices, orthodontic splints

Tips and tricks:

Vacuum-sealed material, does not require predrying. PMMA adjustment can be done without bonding agent, use of ISOFOLAN® required.

Available in:

1.2 / 1.8 / 2.5 / 3.0 / 4.0 mm
x 125 mm Ø

Technical data

Material:	TPU
Hardness:	35 Sh. D
Tensile strength:	39 MPa
Impact strength:	no fracture
Material:	PETG
E-modulus:	2,200 MPa
Hardness:	78 Sh. D
Tensile strength:	53 MPa
Impact strength:	10 KJ/m ²

CE 0044

DURASOFT® blue pd

Indication

Hard-elastic, blue-transparent splint with soft inner surface in particular for easy production of coloured intra-oral anti-snoring appliances. Blue colour increases aesthetics.

Tips and tricks:

Vacuum-sealed material, does not require predrying, PMMA adjustment can be done without bonding agent, use of ISOFOLAN® foil required.

Available in:

1.8 / 2.5 mm x 125 mm Ø

Technical data

Material:	TPU
Hardness:	35 Sh. D
Tensile strength:	39 MPa
Impact strength:	no fracture
Material:	PETG
E-modulus:	2,200 MPa
Hardness:	78 Sh. D
Tensile strength:	53 MPa
Impact strength:	10 KJ/m ²

CE

DURASOFT® seal

Indication

Discolouration resistant inner layer for DURASOFT® pd splints. Specifically designed for long-term use splints such as intra-oral anti-snoring devices.

Tips and tricks:

Can be combined in the lamination process with DURASOFT pd. Use of ISOFOLAN® foil required.

Available in:

0.1 mm x 125 mm Ø
clear, blue, magenta

Technical data

Material:	Thermoplastic polyurethane
E-modulus:	1,900 MPa
Hardness:	82 Sh. D
Tensile strength:	69 MPa
Impact strength:	no fracture

CE 0044

HARDCAST®

Indication

Copings, spacer foil for bleaching splints

Tips and tricks:

Copings can be thermoformed with spacer foil.

Available in:

0.4 / 0.6 / 0.8 mm x 125 mm Ø

Technical data

Material:	Polypropylen
E-modulus:	890 MPa
Hardness:	Rockwell 72
Tensile strength:	28 MPa
Impact strength:	600 KJ/m ² o. no fracture

IMPRELON® opaque

Indication

Individual impression trays

Tips and tricks:

Ready-made handles can be attached using resin, use moistened foam discs as spacer, edges can be smoothed with monomer.

Available in:

2.0 / 3.0 mm x 125 mm Ø
opaque

Technical data

Material:	Polystyrol
E-modulus:	2,500 MPa
Tensile strength:	18 MPa

IMPRELON® clear

Indication

Individual impression trays

Tips and tricks:

Ready-made handles can be attached using resin, use moistened foam discs as spacer, edges can be smoothed with monomer.

Available in:

2.0 / 3.0 mm x 125 mm Ø
clear

Technical data

Material:	Polystyrol
E-modulus:	3,200 MPa
Hardness:	Rockwell 150
Tensile strength:	55 MPa
Impact strength:	14 KJ/m ²

IMPRELON® white

Indication

Bite registration tray

Tips and tricks:

Ready-made handles can be attached using resin, use moistened foam discs as spacer, edges can be smoothed with monomer.

Available in:

2.0 / 3.0 mm x 125 mm Ø
white

Technical data

Material:	Polystyrene
E-modulus:	1,850 MPa
Hardness:	Rockwell 80
Tensile strength:	17 MPa
Impact strength:	> 30 KJ/m ²

IMPRELON® S pd

Indication

Miniplast splints, retention splint,
OSAMU-Retainer®

Tips and tricks:

Can be used in conjunction with the bonding agent OSAMU-Bond, BIOPLAST® 1.5 mm and IMPRELON® S pd 0.75 mm for the production of an OSAMU-Retainer®.

Available in:

0.5 / 0.75 / 1.0 / 1.5 / 2.0 / 3.0 mm
x 125 mm Ø

Technical data

Material:	Copolyester
E-modulus:	1,624 MPa
Hardness:	78 Sh. D
Tensile strength:	49 MPa
Impact strength:	no fracture

ISOFOLAN®

Indication

Best suited as insulation foil when pressurising hardelastical material for splints.

Tips and tricks:

The ISOFOLAN® foil should be thermoformed, cut and perforated interdentially prior to the fabrication of any kind of hard splints.

Available in:

0.1 mm x 125 mm Ø
natur

Technical data

Material:	HD PE
E-modulus:	600 MPa
Hardness:	45 Sh. D
Tensile strength:	30 MPa
Impact strength:	no fracture

COATING FOIL

Indication

Packaging or protection foil for model dispatch or study models.

Tips and tricks:

Can be used together with the Skin packaging card board for an attractive packaging.

Available in:

0.15 mm x 125 mm Ø

Technical data

Material:	Poly-vinyl-chloride
Hardness:	82 Sh. D
Tensile strength:	46 MPa
Impact strength:	> 550 KJ/m ²

SPACER FOIL

Indication

Space maintainer and insulation foil for COPYPLAST® and HARDCAST®, for the fabrication of copings made of HARDCAST®.

Tips and tricks:


Can be used with HARDCAST® or COPYPLAST® for delicate dies.

Available in:

0.1 mm x 125 mm Ø
clear

Technical data

Material:	Poly-vinyl-chloride
Hardness:	78 Sh. D
Tensile strength:	42 MPa
Impact strength:	600 KJ/m ²

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SCHEU-DENTAL GmbH