

Cleansing and Care

for orthodontic appliances, splints and dentures



Dissolve the content of a bag with 150 ml of water.



Put the orthodontic appliance/splint into the cleansing agent and leave for 30 minutes.



Any remaining sediment can be removed with toothpaste and toothbrush.



The care spray refreshes the appliance before insertion.

Due to contact with food, drinks, pharmaceuticals and individual mouth environment, splints, orthodontic appliances and dentures may change their colour, taste and smell. Sediments can stick to the surface. With CETRON® you or your patient purchase a good cleansing agent with a fast cleansing and long lasting care effect. The good value cleansing solution extends the time of comfortable use and leads to a high patient compliance. All denture and orthodontic devices can be cleaned, freed from sediments and refreshed.

Delivery Program

// CETRON® Care Set	1 box with 5 bags (15 g)	3236.1
// Refill package	5 bags (15 g)	3239.1
// Refill package	25 bags (15 g)	3239.2
// Refill package	250 bags (15 g)	3239.3
// CETRON® Care spray	75 ml	3238.1

