



With confidence an attractive smile

At first glance, CA® CLEAR ALIGNER is just another aligner: Made from clear-transparent material, almost invisible, pleasantly smooth and comfortable. At second glance, many things are different with CA® CLEAR ALIGNER: Different material thicknesses and an intelligent splint design for faster treatment results, the optional insertion of a specialised expansion screw, complete treatment control thanks to step-by-step planning.

CA® CLEAR ALIGNER – one system, many options.

CA CLEAR ALIGNER

Contact

SCHEU-DENTAL GmbH
Am Burgberg 20
58642 Iserlohn
GERMANY

Phone: +49 23 74 92 88-0
Fax: +49 23 74 92 88-90
E-Mail: service@scheu-dental.com
www.ca-clear-aligner.com

CA DIGITAL GmbH
Willettstraße 10
40822 Mettmann
GERMANY

Phone: +49 21 04 83 37 12-0
Fax: +49 21 04 83 37 12-3
E-Mail: info@ca-digit.com
www.ca-digit.com

GB 5.000/05/16 G REF PM006402



CA CLEAR ALIGNER



A beautiful smile at any time.



CA DIGITAL

CA DIGITAL is your direct contact and partner to get acquainted with the digital splint therapy. We'll assist you with all questions concerning precise, digital treatment planning, offering you individual cooperation options.

CA DIGITAL delivers you printed 3D set-up models for your splint fabrication at site. This means control of treatment is completely in your hands.

Detailed examples for the possible indications can be found on www.ca-clear-aligner.com. Here you can download our newsletters showing many cases selected by the clinical board of advisers as well as new patient cases and developments in the treatment concept presented by Dr. Pablo Echarrri at a regular base.



The CA® CLEAR ALIGNER therapy concept is versatile: Space closure, crowding, deep bite, cross bite ... these are only a few examples for its effectiveness when it comes to various kinds of tooth corrections.

Individual and precise correction of malocclusion.



CA LAB


Choosing the CA LAB option, all treatment steps from planning, set-up fabrication until realisation of the splint will be carried out in your practice resp. your laboratory.

The necessary equipment consists of a pressure moulding machine type MINISTAR S® or BIOSTAR®, CA® material in soft, medium and hard as well as the CA® SMART software including webcam, CA® photo support and CA® CHECKER.

The convincing aligner therapy made by you.

The proven and successful CA® CLEAR ALIGNER system has been developed by SCHEU-DENTAL especially for treating adults and allows aligning teeth easily and aesthetically.

More than 100.000 patients have been successfully treated with CA® CLEAR ALIGNER these last 4 years only – this means valuable experience resulting from this impressive number. Our expert committee consisting of clinicians and ortho technicians is continuously exchanging this experience to the benefit of CA® users worldwide.



■ Countries with certified CA® partners

The right concept for each practice.

Many things are different with CA® CLEAR ALIGNER, e.g. the form of the splints or the material thicknesses. But what makes CA® CLEAR ALIGNER really different: It's up to you to decide. Do you have to react to modifications during treatment? Do you want to insert additionally a VECTOR® 40 expansion screw in order to correct crowdings? And above all: Do you want to carry out the complete fabrication process on your own, to cooperate with a certified laboratory or to have the treatment planned and realised by cutting-edge 3D technology? CA® CLEAR ALIGNER offers you the choice.




Added value is completely in your hands.

Innovative 3D technology for the most precise treatment results.

The 3-Splint-System.

Three material thicknesses, worn during four weeks for a unique and beautiful smile: With each treatment step, the patient will be wearing three splints of different thicknesses. Thanks to the steady increase in pressure, teeth are aligned efficiently, yet smoothly and with extraordinary comfort. The special splint design with additional covering of the gingival margin (2-3 mm) increases power transmission, providing even more efficiency in treatment.



Your advantages:

- Added value is completely in your hands
- Complete treatment control
- Treatment corrections are possible at any time
- Intelligent software CA® SMART, digital measurement technique for set-up, visual data transfer in real time by webcam

- Option to cooperate with certified CA® laboratories or CA DIGITAL
- Practical training courses in the SCHEU-ACADEMY

CA LAB

Advantages for your patient:

- Invisible malocclusion correction
- Simple care, can be taken out for cleansing
- Transparent splint without brackets or wires
- High comfort, pleasant feeling in the mouth
- High quality, biocompatible material
- Speaking and laughing without problems
- Short treatment time
- Clinical experience



CA DIGITAL

Your advantages:

- Precise, digital 3D technology for planning and set-up fabrication
- Teeth correction based on the clinical CA® records by Dr. Pabo Echarri
- Step-by-step procedure or up to 3 steps per impression
- Treatment corrections are possible at any time



- Unique patient models, made by 3D printing in anatomical colour
- Situation and final models in anatomical colour help to motivate your patients for the CA® CLEAR ALIGNER treatment
- Individual cooperation: It's up to you – choose full or only partial services, STL data or models