



**BIOPLAST®**  
**XTREME**

## Innovative Construction

First pressure-moulding material with two different degrees of hardness

### 1 Most comfortable

This custom mouth guard has the perfect fit. During sports, the slight impression worked into the mouth guard stabilises the teeth in the lower jaw.



**BIOPLAST®**  
**XTREME**

**BIOPLAST®**  
**XTREME PRO**

### 2 Can be combined

Your braces treatment isn't affected, as brackets can be blocked out.

### 3 Highly effective

Wearing a mouth guard significantly reduces the risk of injury to lips, teeth, jaws, and temporomandibular joint.

### 4 Impact absorbing

The impact is absorbed elastically and distributed with less force over a time delay.

### 5 Hygienic

**CETRON®** cleansing powder, toothbrush and toothpaste are suitable for gentle care.

## Precision-Fit Mouth Guard

in just a few steps.

### 1 Consult with your dentist

Depending on the types of sports and recreational activities practiced, the dentist will advise you or your child on selecting the optimal mouth guard.

### 2 Impression

An impression of the upper jaw is the basis for sports mouth guards.

### 3 Producing the mouth guard

It is custom-made with precision in the laboratory. It covers the upper row of teeth and has a slight impression for the lower teeth.

### 4 Fitting and examination

The fit will be tested regularly by your dentist and adjusted as needed on an individual basis.

# Perfection is protecting a smile.



The safe, ideal mouth guard for any application.

# BIOPLAST®

Provided by:

[www.scheu-dental.com](http://www.scheu-dental.com)

**SCHU**

GB 04/2.000/01/15 G REF 0308.02



**SCHU**  
Dental Technology



# BIOPLAST®

What is it?

BIOPLAST® is a high-quality ethylene vinyl acetate (EVA) distinguished by its first-class impact absorption and outstanding elasticity despite its extremely light weight.

BIOPLAST® mouth guards for active sports are custom-made just for you by dental technicians, who work together with your dentist. This guarantees that your personal mouth guard will suit your particular needs for all types of sports.



## For any challenge

Three mouth guard types




### BIOPLAST® MASTER

- proven, two-layer mouth guard
- for all recreational sports
- double-layer laminate construction for the best protection and comfort
- 10 colour variations are available



### BIOPLAST® XTREME




### BIOPLAST® XTREME

- with integrated reinforcement section in the front teeth area
- suited for contact sports, as it protects against hard and concentrated hits
- double-layer laminate construction for the best protection and comfort



### BIOPLAST® XTREME PRO



### BIOPLAST® XTREME PRO

- with integrated reinforcement area and supplemental hard-elastic insert around the front teeth area.
- three-layer material structure for maximum protection and comfort
- for professional contact sports with expected concentrated, hard hits (impact is distributed across a greater area)



“A large number of sports accidents include injuries to the teeth and the periodontium.”