

<b>REF 3370 HM Tricutter fine (green)</b> Finishing of hard elastic foils and splints	<b>Max. 7.500 r/min</b>
<b>REF 3381 DIMO®PRO</b> Final trimming and polishing	<b>Max. 10.000 r/min</b>
<b>REF 5262 CA®-Carbide Bur</b> Trimming of plaster segments	<b>Max. 20.000 r/min</b>
<b>REF 5288 Diamond separating disc</b> Cutting of plaster segments with a smooth cut	<b>Max. 5.000 r/min</b>
<b>REF 5293 Set-Up separating disc</b> Ø 40 mm, 0.5 mm Cutting of plaster segments with a rough cut	<b>Max. 15.000 r/min</b>
<b>REF 5371 Set-Up model trimmer</b>	<b>Max. 10.000 r/min</b>