

3
different
versions

The Bite Fix® bite registrations


The Bite Fix® bite registrations enable simple and easy construction bites in snoring therapy. The decentralized handle ensures an easy midline control and comes with fixing screw and mm-scale for an uncomplicated control of the protrusion adjustment.

Three versions are available:


Bite Fix® – for multiple use

Bite Fix® Scan – for any model scanner


Bite Fix® One – for single use



Bite Fix®



Bite Fix® Scan



Bite Fix® One

These are the common features of our Bite Fix® bite registrations:

- // Decentralized handle guarantees an easy midline control.
- // Fixing screw ensures smooth handling.
- // mm scale enables an easy protrusion adjustment.
- // The material is torsion-resistant and light, yet durable.
- // The delicate design allows for a handy operation.
- // Thanks to the fixed vertical opening, TAP® and IST devices can be easily realised.

Die Bite Fix®

bite registrations

Bite Fix® – Quick. Safe. Precise.

Combines all the advantages:

- // mm-scale to control the protrusion adjustment
- // Fixing screw for smooth handling
- // Decentralised handle for an easy midline control
- // Torsion-resistant

Plus: The Bite Fix® bite registration is made of sterilisable plastic, the basic element (red) is for multiple use.



Bite Fix® kit



Bite Fix® bite fork

Delivery program

// Bite Fix® bite registration	1 kit	5474
// Bite Fix® bite fork	10 pcs.	5475

Bite Fix® Scan – Scanning models in protrusion.

Offers even more:

- // mm-scale to control the protrusion adjustment
- // Fixing screw for smooth handling
- // Decentralised handle for an easy midline control
- // Basic element (red) made of sterilisable plastic
- // Torsion-resistant

Plus: Thanks to the pre-determined breaking point on the bite fork, the Bite Fix® Scan fits into any model scanner.



Bite Fix® Scan kit



Bite Fix® Scan bite fork

Delivery program

// Bite Fix® Scan bite registration	1 kit	5492
// Bite Fix® Scan bite fork	10 pcs.	5493

Bite Fix® One – For single use.

With focus on the essential:

- // mm-scale to control the protrusion adjustment
- // Fixing screw for smooth handling
- // Decentralised handle for an easy midline control
- // Torsion-resistant

Plus: Although less costly, the Bite Fix® One combines all relevant features for an easy handling.



Bite Fix® One kit

Delivery program

// Bite Fix® One bite registration	1 kit	5503
------------------------------------	-------	------

CE

